

EpiFaith® products are proudly recognized by major healthcare purchasing networks and institutions

VA



U.S. Department  
of Veterans Affairs


vizient.

PREMIER



# Safer Epidural Starts Here.

**EpiFaith® Smart Syringe**



**Smart Signal**  
(After LOR)



Scan to learn more.

# SAFE. SIMPLE. SMARTER EPIDURALS.




## See the Pressure, Enhance Safety

EpiFaith® uses an automated pressure-sensing system that delivers real-time and continuous detection—providing an objective, visual signal upon entry into the epidural space, as validated by global clinical use for enhanced safety.

## Teaching Made Easy

Clear, visual feedback simplifies training and improves communication. EpiFaith helps learners focus on what matters—precise needle placement.



All connector types—  
Luer Slip, Luer Lock, and **NRFit®**—are available

## Control Made Simple

With two-handed control and real-time pressure visuals, EpiFaith empowers providers to place needles with greater stability, accuracy, and safety.



## ENDPOINT UNDER EVALUATION

### wet tap rate ↓

*Retrospective and presentation data reported lower wet tap incidence.*

Poster, TSA 2024.  
BWH presentation, SOAP 2025.

### assistance needed ↓

*Early comparative observations showed a downward trend among junior clinicians.*

Brigham and Women's Hospital data presented at ASA 2025.

### false positive ↓

*A published study noted fewer false-positive LOR events.*

Pain physician, 28(4), E431–E435.

*For more information on clinical results, please contact your local sales representative.*

USA			CE		
Type	Ref. #	Qty/Case	Type	Ref. #	Qty/Case
Luer Slip	FM-01SLR	50	Luer Slip	EF-01SLR	50
Luer Lock	FM-01SLL		Luer Lock	EF-01SLL	
NRFit	FM-025NR		NRFit	EF-025NR	



## Taiwan

(+886) 2-2567-2959  
12F-2, No. 206, Sec. 2, Nanjing E. Rd.,  
Zhongshan Dist., Taipei 104, Taiwan

## US

(+1) 737-237-8337  
5900 Balcones Drive Ste 25706  
Austin, TX 78731, USA

## EU

(+31) 33 445 12 39  
Chromiumweg 139, 3812 NM, Amersfoort,  
Utrecht, Netherlands