

EpiFaith® products are proudly recognized by major healthcare purchasing networks and institutions

VA



U.S. Department
of Veterans Affairs


vizient.

PREMIER



Safer Epidural Starts Here.

EpiFaith® Smart Syringe



Smart Signal
(After LOR)



Scan to learn more.

SAFE. SIMPLE. SMARTER EPIDURALS.




See the Pressure, Enhance Safety

EpiFaith® uses an automated pressure-sensing system that delivers real-time and continuous detection—providing an objective, visual signal upon entry into the epidural space, as validated by global clinical use for enhanced safety.

Teaching Made Easy

Clear, visual feedback simplifies training and improves communication. EpiFaith helps learners focus on what matters—precise needle placement.



All connector types—
Luer Slip, Luer Lock, and **NRFit®**—are available

Control Made Simple

With two-handed control and real-time pressure visuals, EpiFaith empowers providers to place needles with greater stability, accuracy, and safety.



PROVEN REDUCTIONS

wet tap rate

90.9% ↓

*compared to
traditional devices*

assistance
needed

31.8% ↓

*by newly qualified
physician, compared
to traditional devices*

false
positive

60.0% ↓

*with saline, compared
to literature reports*

For more information on clinical results, please contact your local sales representative.

USA			CE		
Type	Ref. #	Qty/Case	Type	Ref. #	Qty/Case
Luer Slip	FM-01SLR	50	Luer Slip	EF-01SLR	50
Luer Lock	FM-01SLL		Luer Lock	EF-01SLL	
NRFit	FM-02SNR		NRFit	EF-02SNR	

Taiwan

(+886) 2-2567-2959
12F-2, No. 206, Sec. 2, Nanjing E. Rd.,
Zhongshan Dist., Taipei 104, Taiwan

US

(+1) 737-237-8337
5900 Balcones Drive Ste 25706
Austin, TX 78731, USA

EU

(+31) 33 445 12 39
Chromiumweg 139, 3812 NM, Amersfoort,
Utrecht, Netherlands